


## Catalog Number

## 58820-2L-CS

**Bottles, Mobile Phase Reservoir, safety, coated glass, UV cut off of 280nm, clear, 2 L, screw Top, GL45 threads. Caps sold separately. ARE-Applied Research brand. CLEARANCE. 5/CS.**



### Additional Info

This Clear, Plastic Coated, Borosilicate Glass Bottle has a 2L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir.

These bottles are in excellent condition and are not "seconds", damaged, returned goods or defective in any way.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

## \$300.94

Printed On 2024/04/30