

## Catalog Number

58905-S-K

**Syringe, Chrom. Disposable, 5ml, Luer Lock, Non Sterile.**  
MICROSOLV Brand. Bulk Packed. 1000/CS.



### Additional Info

These Luer Lock Tip, *Chrom Syringes* are chosen for many uses in chromatography as well as used with syringe filters or for any lab situation needing an inert, non-reactive, plastic syringe. The Luer Lock device provides additional safety by securing the syringe with the needle or filter to prevent the devices from popping off the syringe during times of higher pressures.

Latex free, they do not contain rubber, silicone oil, styrene or DEHP that is found in other plastic syringes. These all plastic syringes are manufactured from laboratory grade polypropylene and polyethylene. More chemically resistant than syringes with

**Price**  
**\$527.43**

## Additional Info

rubber-tipped plungers, they are also equipped with a positive "click", safety stop to help prevent accidental spills. Bulk packed and not individually wrapped, saves time unwrapping and waste disposal costs.

Click [HERE](#) for smaller package quantity or click [HERE](#) for standard package quantity.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$527.43**