

Catalog Number

59001-25P

HPLC Column, MicroBee Silica, 10um, 4.6mm ID x 250mm Length, 125A. Cogent Brand. 1 EA.



Additional Info

This Very Popular Cogent MicroBee Silica™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is irregular shaped silica that is not modified in any way. This makes it an ideal Alternate in older USP Methods that use the L3 Classification for columns.

Often used for hydrophilic or polar compounds in Normal Phase HPLC.

Price
\$781.79

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$781.79