

Catalog Number

59052-02-4

HPLC Column, MicroBee Phenyl, 10um, 4.6mm ID x 300mm Length, 125A. Cogent Brand. 1 EA.



Additional Info

This Cogent MicroBee Phenyl™ Stationary Phase is part of the Similars™ product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A", 10um, 125A Stationary Phase is fully end capped and is manufactured with a 4 carbon ligand that has an attached Phenyl group. Bonded to irregular shaped silica, this makes it an ideal Alternate in USP Methods that use the L11 Classification and Microbondapak™ or similar Phenyl Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$764.59

Printed On 2024/04/29