

## Catalog Number

59150-01

**HPLC Column, MicroBee C18**, 10um, 3.9mm ID x 150mm Length,  
125A. Cogent Brand. 1 EA.



### Additional Info

This Cogent MicroBee C18™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to irregular shaped silica. This makes this column an ideal Alternate in USP Methods that use the L1 Classification and Microbondapak™ or similar C18 columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

*This is one of our Most Popular MicroBee C18 Columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$691.26**

Printed On 2024/07/27