

## Catalog Number

59300-03

**HPLC Column, MicroBee Cyano, 10um, 3.9mm ID x 300mm Length, 125A. Cogent Brand. 1 EA.**  
MOST POPULAR.



### Additional Info

This Cogent MicroBee Cyano™ Stationary Phase is part of the Similar™ product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A", 10um, 125A Stationary Phase is not end capped and is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond. Bonded to irregular shaped silica, this makes it an ideal Alternate in USP Methods that use the L10 Classification and Microbondapak™ or similar Cyano Columns.

**Price**

**\$764.59**

## Additional Info

*This is one of our Most Popular MicroBee Cyano Columns.*

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### **Price**

**\$764.59**