

## Catalog Number

59300-04

**HPLC Column, MicroBee Amino**, 10um, 3.9mm ID x 300mm Length,  
125A. Cogent Brand. 1 EA.  
MOST POPULAR.



### Additional Info

This Cogent MicroBee Amino™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is not end capped and manufactured with traditional bonding of a 3 carbon chain ligand with a terminal Amino group, to irregular shaped silica. This Column is an ideal Alternate in USP Methods that use the L8 Classification and Microbondapak™ or similar Amino columns.

*This Column is packed in and shipped in Normal Phase Solvents. This is one of our Most Popular MicroBee Amino Column.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$764.59**

Printed On 2024/05/08