

Catalog Number

62508-25P

HPLC Column, Hyper Similar C8 Octyl, 5um, 4.6mm ID x 250mm
Length, 150A. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent Hyper Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Thermo Fisher Hypersil™ Columns or Similar C8 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids but expect less retention than with a C18 Column.

This is our Most Popular Hyper Similar C8™ Column.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/04/30