

# Catalog Number

62518-10P-3

**HPLC Column, Hyper Similar C18, 5um, 3.0mm ID x 100mm Length, 150A. Cogent Brand. 1 EA. CUSTOM ORDER.**



## Additional Info

This Cogent Hyper Similar C18™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years. The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L1 Classification and Thermo Fisher Hypersil™ Columns or Similar C18 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl

**Price****\$785.53**

## Additional Info

Phase is excellent for long chain fatty acids.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$785.53**