

Catalog Number

63098-10P-W

HPLC Column, C18 ODS 3, 5um, 4.6mm ID x 100mm Length, 100A with WCS hardware. Partisil Brand. 1 EA.



Additional Info

This Partisil C18™ Octadecyl Silane (ODS3), Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is fully end capped and modified with traditional bonding of a C18 ligand. This Column is written in many USP Methods; Partisil Silica used in older methods with the L1 Classification.

Price**\$1,445.15**

Additional Info

Often used for hydrophobic or non polar compounds in Reverse Phase HPLC, this column is supplied with the Whatman Compression Sealing System RAC-II (WCS).

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,445.15