

## Catalog Number

63118-25P-W

**HPLC Column, C18 ODS 1**, 10um, 4.6mm ID x 250mm Length, 100A with WCS hardware. Partisil Brand. 1 EA.



### Additional Info

This Partisil C18™ Octadecyl Silane (ODS 1), Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is not end capped but modified with traditional bonding of a C18 ligand. This Column is written in many USP Methods; Partisil Silica used in older methods with the L1 Classification.

*Often used for hydrophobic or non polar compounds in Reversed Phase HPLC, this column is supplied with the Whatman Compression Sealing System RAC-II (WCS).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,533.76**

Printed On 2024/05/05