

## Catalog Number

63185-C2

**HPLC Column, SAX**, 10um, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisil Brand. 1 EA.



### Additional Info

This Partisil SAX™ Strong Anion Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is modified with a ligand that contains tetra-methyl ammonium. This Column is written in many USP Methods that call for a Partisil SAX Column and the L23 Classification.

*Often used for purification or analysis of biological compounds in Ion Exchange HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires end fitting kit [66900-WEF](#).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,399.75**

Printed On 2024/05/03