

Catalog Number

63189-25P-W

HPLC Column, SCX (S-P), 10um, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisil Brand. 1 EA.



Additional Info

This Partisil SCX™ (S-P) Strong Cation Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is modified with a ligand that contains sulfonic acid. This Column is written in many USP Methods that call for a Partisil SCX Column and the L22 Classification.

Often used for purification or analysis of biological compounds in Ion Exchange HPLC, this column is supplied with the Whatman Compression Sealing System (WCS).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,549.25

Printed On 2024/05/03