

# Catalog Number

63189-25P-W

**HPLC Column, SCX (S-P)**, 10um, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisil Brand. 1 EA.



## Additional Info

This Partisil SCX™ (S-P) Strong Cation Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is modified with a ligand that contains sulfonic acid. This Column is written in many USP Methods that call for a Partisil SCX Column and the L22 Classification.

**Price****\$1,549.25**

## Additional Info

*Often used for purification or analysis of biological compounds in Ion Exchange HPLC, this column is supplied with the Whatman Compression Sealing System (WCS).*

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### **Price**

**\$1,549.25**