

Catalog Number

64418-10P-2

HPLC Column, Nova Similar C18, 4um, 2.1mm ID x 100mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent Nova Similar C18™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Nova-Pak®. This makes this column an ideal Alternate in USP Methods that use the L1 Classification and Waters Nova-Pak™ Columns or Similar C18 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids

A very popular column size, this narrow-bore column provides short run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rate with lower volumes.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$848.36

Printed On 2024/05/05