

Catalog Number

64418-PG2

Guard Column Cartridges, Nova Similar C18, 4um. 4.0mm x 10mm
Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridges are packed with 4µm, Cogent Nova Similar C18™ Stationary Phase. They are end capped (*fully*), with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with 200A Pore Size.

Designed to protect your 3.0, 4.0 or 4.6mm ID, Cogent Nova Similar C18, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$767.01

Printed On 2024/05/03