

## Catalog Number

64424-10P-4

**HPLC Column, Nova Similars Phenyl, 4um, 4.0mm ID x 100mm**  
Length. Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Nova Similars Phenyl™ Stationary Phase is part of the Similars product line that has been a premium value used by HPLC laboratories for many years.

This "Type-B" Stationary Phase is fully end capped along with bonding of a 4 carbon ligand with a terminal Phenyl group to spherical shaped silica that has very Similar surface area to the original Nova-Pak®. This makes this column an ideal Alternate in USP Methods that use the L11 Classification and Waters Nova-Pak™ Columns or Similar Phenyl Columns.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$785.53**

Printed On 2024/05/05