

## Catalog Number

65308-15P-2

**HPLC Column, Zx Similar C8 Octyl**, 3um, 2.1mm ID x 150mm Length.  
Cogent Brand. 1 EA.



### Additional Info

This Cogent Zx Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Agilent's Zorbax® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Agilent's Zorbax™ brand C8 HPLC Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids with less Retention than a C18 Column.

*This narrow-bore column provides good run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS when additional resolution is desired.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,015.95**

Printed On 2024/05/01