

Catalog Number

66000-C1

HPLC Column, Silica, 5um, 4.6mm ID x 125mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.



Additional Info

This Partisphere Silica™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is not modified in any way. This Column is written in many USP Methods; Partisil Silica used in older methods with the L3 Classification.

Often used for hydrophilic or polar compounds in Normal Phase HPLC, this column is

Price**\$1,056.08**

Additional Info

the cartridge style, Whatman Void Sealing System (WVS) and requires end fitting kit [66900-WEF](#) which is sold separately.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,056.08