

# Catalog Number

66018-C1

**HPLC Column, C18, 5um, 4.6mm ID x 125mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.**



## Additional Info

This Partisphere C18™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is bonded with a ligand that is 18 carbons long. This Column is written in many USP Methods; Partisil Silica used in older methods with the L1 Classification.

*Often used for hydrophobic or non polar compounds in Reversed Phase HPLC, this*

**Price****\$1,083.88**

## Additional Info

*column is the cartridge style, Whatman Void Sealing System (WVS) and requires end fitting kit [66900-WEF](#) which is sold separately.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,083.88**