

# Catalog Number

66086-C2

**HPLC Column, SCX, 5um, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.**



## Additional Info

This Partisphere SCX™ Strong Cation Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is modified with a ligand that contains sulfonic acid. This Column is written in many USP Methods that call for a Partisil SCX Column and the L22 Classification.

*Often used for purification or analysis of biological compounds in Ion Exchange*

**Price****\$1,153.37**

## Additional Info

*HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires the end fitting kit [66900-WEF](#) which is sold separately.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,153.37**