

Catalog Number

66087-C2

HPLC Column, PAC, 5um, 4.6mm ID x 250mm Length, 100A with WVS
Cartridge hardware. Partisphere Brand. 1 EA.



Additional Info

This Partisphere PAC™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is not end capped and is modified with a ligand that contains both Amino and Cyano functional groups. This Column is written in many Methods that call for a Partisil Amino-Cyano (PAC) Column.

This column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires the end fitting kit [66900-WEF](#) which is sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,127.99

Printed On 2024/05/02