

Catalog Number

66090-P25-W

HPLC Column, WAX Max-1, 5 μ m, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.



Additional Info

This Partisphere WAX™ Weak Anion Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped, and is modified with a ligand that contains proprietary ionic groups. This Column is written in many USP Methods that call for a Partisphere WAX Column and the L23 Classification.

Often used for purification or analysis of biological compounds in Ion Exchange

Price**\$1,320.07**

Additional Info

HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires the end fitting kit [66900-WEF](#) which is sold separately.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,320.07