

Catalog Number

66900-PAC-GC

PAC Guard Cartridges, for use with WCS or WVS holders with Partisphere HPLC columns, 5/PK. CUSTOM ORDER.



Additional Info

These guard column cartridges are packed with 5um, PAC stationary phase to be used with Partisphere brand HPLC Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$450.41

Printed On 2024/05/05