

Catalog Number

66900-PAC-GC

PAC Guard Cartridges, for use with WCS or WVS holders with Partisphere HPLC columns, 5/PK. CUSTOM ORDER.



Additional Info

These guard column cartridges are packed with 5um, PAC stationary phase to be used with Partisphere brand HPLC Columns.

Price
\$460.55

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$460.55