



# Catalog Number

# 67585-15P-2

**HPLC Column, IEX SAX, 5um, 2.1mm ID x 150mm Length, 100A.**  
Cogent Brand. 1 EA.



## Additional Info

This Cogent Ion Exchange Column is a Strong Anion Exchanger (SAX has a positively charged surface ) Stationary Phase that has been a premium value and used by HPLC laboratories for many years. It is manufactured with Quaternary Amine functional groups that are always charged.

We use only the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L9 Classification or for acids analysis or purification.

**Price**

**\$915.06**



## Additional Info

With a 5 $\mu$ m Particle, this column works well and produces excellent data for challenging IEX applications that include Anions (negatively charge compounds).

### MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$915.06**