

## Catalog Number

67586-15P

**HPLC Column, IEX SCX, 5um, 4.6mm ID x 150mm Length, 100A.**

Cogent Brand. 1 EA.

**MOST POPULAR.**



### Additional Info

This Very Popular Cogent Ion Exchange Column is a Strong Cation Exchanger (SCX has a negatively charged surface) Stationary Phase that has been a premium value and used by HPLC laboratories for many years. It is manufactured with Sulfonic Acid functional groups that are always charged.

We use only the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L9 Classification or for "Bases" analysis or purification.

With a 5µm Particle, this column works well with challenging IEX applications that include Cations (positive charge compounds).

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$785.53**

Printed On 2024/05/18