

Catalog Number

68004-P30-50

Prep HPLC Column, RP C4, 10um, 100A. 30mm ID x 50mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 carbon (*butyl*) ligand to the highest purity silica available.

With a 10 µm Particle, this column works well with higher flow rates and lower back pressures. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for preparative HPLC for small molecule purification. Expect Less retention than with a C18 column but more than with a C8.

Price**\$4,195.88**

Additional Info

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$4,195.88