

## Catalog Number

68004-P50-150

**Process HPLC Column, RP C4**, 10 $\mu$ m, 100A. 50mm ID x 150mm Length. Cogent Brand. 1 EA.



### Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 carbon (*butyl*) ligand to the highest purity silica available.

With a 10  $\mu$ m Particle, this column works well with higher flow rates and lower back pressures. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for in process purification of small molecules. Expect Less retention than with a C18 column but more than with a C8.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$16,382.00**

Printed On 2024/04/25