

Catalog Number

68004-P50-1503M

Process HPLC Column, RP C4, Wide Pore, 10um, 300A. 50mm ID x 150mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available. This makes this column ideal for large molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used in process for purification of proteins, peptides, enzymes and other large, hydrophobic or neutral compounds in Reversed Phase preparative HPLC, this Butyl Phase is excellent for bio molecules.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$16,382.00

Printed On 2024/04/25