

Catalog Number

68004-SP100-3M

Semi Prep HPLC Column, RP C4 Wide Pore, 10um, 300A. 10mm ID x 100mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. The diameter and length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of proteins, peptides, enzymes and other biological compounds making it ideal for Semi Preparative HPLC.

Price**\$1,680.56**

Additional Info

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,680.56