

## Catalog Number

68008-05P-4-3M

**HPLC Column, RP C8 Octyl, Wide Pore**, 10um, 4.0mm ID x 50mm  
Length, 300A. Cogent Brand. 1 EA.  
*CUSTOM ORDER.*



### Additional Info

This Cogent RP C8™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L7 Classification, for large molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for large hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is also excellent for long chain fatty acids.

*This column is often chosen as a scouting column for preparative methods.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$576.05**

Printed On 2024/05/07