

Catalog Number

68008-P21-PG1-3

Prep Guard Column Cartridge, RP C8 Octyl, Wide Pore, 10um, 300A.
21.2mm ID x 10mm Length. Cogent Brand. 1 EA.



Additional Info

This Guard Column Cartridge is packed with 300Å, 10µm Cogent RP C8™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available. It is designed to protect your RP C8 Wide Pore, Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder. Click [HERE](#) for more information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$447.95

Printed On 2024/05/02