

## Catalog Number

68008-SP250-3M

**Semi Prep HPLC Column, RP C8 Octyl, Wide Pore, 10um, 300A.**  
10mm ID x 250mm Length. Cogent Brand. 1 EA.



### Additional Info

This Cogent RP C8™, 300Å, Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L7 Classification, for large molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. The diameter and long length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of large hydrophobic or neutral compounds, making it ideal for Semi Preparative HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,384.79**

Printed On 2024/05/02