

Catalog Number

68022-P50-50

Process HPLC Column, RP PFP, 10um, 100A. 50mm ID x 50mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

With a 10µm Particle, this column works well with higher flow rates producing low back pressures in preparative HPLC making this column ideal for purification of halogenated or polar samples.

A specialized HPLC System is required for this column which is packed and shipped in Reversed Phase Solvents.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$10,248.15

Printed On 2024/05/05