

Catalog Number

68024-15P-4

HPLC Column, RP Phenyl, 10um, 4.0mm ID x 150mm Length, 100A.
Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP Phenyl™ HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon linkage to a terminal Phenyl Ring known as Propyl Phenyl, to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification.

Used to separate Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC, these Phenyl

Price**\$674.70**

Additional Info

Columns are an excellent choice.

With a 10 μ m Particle, it works well with high flow rates and and produces lower back pressures.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$674.70