

## Catalog Number

68024-P21-PG1

**Prep Guard Column Cartridge, RP Phenyl, 10µm, 100A. 21.2mm ID x 10mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Guard Column Cartridge is packed with 10µm, Cogent RP Phenyl™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available. It is designed to protect your RP Phenyl Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder, *click [HERE](#) for more information about the holder.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/04