

## Catalog Number

68024-P30-100

**Prep HPLC Column, RP Phenyl, 10um, 100A. 30mm ID x 100mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 10 µm Particle works well with higher flow rates and produces lower back pressures than smaller particles.

The extra large diameter and long length of this column provides high loading

## Additional Info

capacity which is great for purification of isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl Rings, making it ideal for Preparative HPLC.

**Notes** - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)