

Catalog Number

68024-P30-50-3M

Prep HPLC Column, RP Phenyl, Wide Pore, 10um, 300A. 30mm ID x 50mm Length. Cogent Brand. 1 EA.



Additional Info

The Cogent RP Phenyl™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories and pilot plants for many years. It is double end capped (fully), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 10 µm particle with 300A pore size works well with higher flow rates and produces lower back pressures than smaller particles and pore sizes.

The extra large diameter and length of this column provides high loading capacity

Additional Info

which is great for purification of larger compounds such as isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Preparative HPLC.

Notes - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com