

## Catalog Number

68024-P50-50

**Process HPLC Column, RP Phenyl, 10um, 100A. 50mm ID x 50mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 10 µm particle works well with higher flow rates and produces lower back pressures than smaller particles.

The extra large diameter and length of this column provides high loading capacity which is great for purification of isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Process HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$9,002.05**

Printed On 2024/05/07