

Catalog Number

68034-05P-3M

HPLC Column, RP Cyano, Wide Pore, 10um, 4.6mm ID x 50mm Length, 300A. Cogent Brand. 1 EA.



### Additional Info

This Cogent RP Cyano™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond and bonded to the highest purity silica available making this Column ideal as an Alternate in USP Methods that use the L10 Classification with larger molecules.

With a 10 µm Particle and 300A Pore Size, this Column works well with higher flow rates and lower back pressures. It is often used for large, polar or neutral compounds

Price

**\$589.01**

## Additional Info

including proteins, peptides, enzymes & natural products in Reversed Phase HPLC. Expect shorter retention times and faster equilibration (*between runs*) than with a non-modified, Silica Column.

*Packed and shipped in Reversed Phase Solvents; click [HERE](#) for Cyano Columns packed in Normal Phase (non polar) Solvents. This Column is used to scout for preparative methods that involve large molecules.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$589.01**