

Catalog Number

68034-25P-3M

HPLC Column, RP Cyano, Wide Pore, 10um, 4.6mm ID x 250mm Length, 300A. Cogent Brand. 1 EA.



Additional Info

This Cogent RP Cyano™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond and bonded to the highest purity silica available making this Column ideal as an Alternate in USP Methods that use the L10 Classification with larger molecules.

With a 10 µm Particle and 300A Pore Size, this Column works well with higher flow rates and lower back pressures. It is often used for large, polar or neutral compounds

Price**\$733.16**

Additional Info

including proteins, peptides, enzymes & natural products in Reversed Phase HPLC. Expect shorter retention times and faster equilibration (*between runs*) than with a non-modified, Silica Column.

Packed and shipped in Reversed Phase Solvents; click [HERE](#) for Cyano Columns packed in Normal Phase (non polar) Solvents.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$733.16