

Catalog Number

68040-05P-3M

HPLC Column, NP Amino, Wide Pore, 10um, 4.6mm ID x 50mm Length, 300A. Cogent Brand. 1 EA.



Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class, wide pore silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods with larger molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups in

Price
\$576.05

Additional Info

Normal Phase or *HILIC* modes, these columns are packed and shipped in Normal Phase Solvents.

This Column is often used for scouting in Preparative Chromatography.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$576.05