

# Catalog Number

68040-P21-250

**Prep HPLC Column, NP Amino, Wide Pore, 10um, 300A. 21.2mm ID x 250mm Length. Cogent Brand. 1 EA.**



## Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 10µm Particle, this Column works well with higher flow rates producing low back pressures compared with smaller particles. Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in

**Price****\$4,450.18**

## Additional Info

Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

**Notes** - *A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$4,450.18**