

Catalog Number

68040-SP150

Semi Prep HPLC Column, NP Amino, 10um, 100A. 10mm ID x 150mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 10µm Particle, this Column works well with higher flow rates producing low back pressures compared with smaller particles. Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,713.04

Printed On 2024/05/07