

Catalog Number

68112-10P

UHPLC Column, RP C12, 1.8 μ m, 100A. 4.6mm ID x 100mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP C12™ UHPLC Stationary Phase and Column has been a premium value and used by HPLC laboratories for several years. It is double end capped (*fully*), and bonded with a rigid Adamantine shaped (diamond like) carbon silane to the highest purity silica available; making this column excellent for high throughput screening (HTS) with rapid equilibration between runs.

The small, 1.8 μ m Particle Size will provide high efficiency even at high flow rates with greater Peak Capacity than other typical Particle Sizes but higher back pressures are normal and to be expected. This column is compatible with any Reversed Phase

Price**\$1,445.78**

Additional Info

mobile phase including 100% water without any loss of retention run to run, which is why this column is often used for high throughput methods that separate complex mixtures. Normally, you can expect less retention than a C18 column but more than a C8 column for the same hydrophobic compounds.

Note: A specialized UHPLC System is required for this high pressure rated column and your samples should be filtered with a 0.2µm Syringe Filter before injection.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,445.78