

## Catalog Number

68304-05P-4

**HPLC Column, RP C4 Butyl, 3 $\mu$ m, 4.0mm ID x 50mm Length, 100A.**  
Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*butyl*) carbon ligand to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L26 Classification.

With a 3 $\mu$ m Particle, this column is more efficient than a 5 $\mu$ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules.

*Expect Less retention than with a C8 column.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$775.06**

Printed On 2024/04/25