

## Catalog Number

68304-10P

**HPLC Column, RP C4 Butyl, 3 $\mu$ m, 4.6mm ID x 100mm Length, 100A.**  
Cogent Brand. 1 EA.  
MOST POPULAR.



### Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*butyl*) carbon ligand to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L26 Classification.

With a 3 $\mu$ m Particle, this column is more efficient than a 5 $\mu$ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$848.36**

Printed On 2024/08/06

## Additional Info

*Expect Less retention than with a C8 column. This is one of our Most Popular C4 columns.*

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$848.36**

Printed On 2024/08/06