

Catalog Number

68324-05P-4-3M

HPLC Column, RP Phenyl, Wide Pore, 3 μ m, 4.0mm ID x 50mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP Phenyl™ Wide Pore HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon chain with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification for large molecules or in new method development.

The 3 μ m Particle makes the column more efficient than a similar column with a 5 μ m Particle, thus you will experience sharper peaks. However, due to the smaller particle size, there can be an increased back pressure so methods that use longer columns

Price

\$792.50

Additional Info

and viscous solvents may exceed the pressure limits of your instrument.

Used to separate large Isomers (*positional*), Natural Products and Related Compounds and some Peptides in Reversed Phase HPLC.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$792.50