

## Catalog Number

68340-10P-2-3M

**HPLC Column, NP Amino, Wide Pore, 3 $\mu$ m, 2.1mm ID x 100mm Length, 300A. Cogent Brand. 1 EA.**



### Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group, to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L8 Classification, for large molecules.

With a 3  $\mu$ m Particle and 300A Pore Size, this column performs well for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups using Normal Phase or *HILIC* modes.

*These columns are packed and shipped in Normal Phase Solvents. This narrow-bore column provides good run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is excellent, and often used with LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$921.70**

Printed On 2024/05/05